




MEMS Transient Analysis by DHM[®]

Transient characterization by Digital Holography Microscope (DHM[®]) encompasses:

- **Excitation signal:** User-programmable waveform shapes
- **Time Domain:** Mega-data-points resolution 4D topography time-sequences
- **Frequency Domain:** Mega-pixels vibration maps, bandwidth up to 200 MHz.
- **Environmental control:** Temperature, pressure, humidity, gas, liquid, etc.

Applications: cMUT, pMUT, inertial sensors, optical and mechanical switches, MEMS microphones, etc.

Time Domain: 4D Mega-resolution Topography



cMUT response to a double period 5 MHz sine burst: Left: Topography color map time-sequence animation. Upper right: 3D perspective. Bottom right: Animated profile cut. (Courtesy to Philips)

DHM[®]'s unique [non-scanning technique](#) measures **4D**, i.e. time-sequences of the successive 3D topographies of a MEMS responding to an excitation signal.

It provides at each time-point:


- **3D surface topography** with a million of data-point resolution
- **profile cut**, giving precise information of surface deformation
- **absolute vertical position** at any point of the measured surface (no need of pre-defined measurement grid).

[Learn more about it](#)

Transients and damping characterization

Lyncée's [stroboscopic module](#) includes a signal generator allowing users to define any waveform shape, in particular **pulses**, **bursts**, and **transient oscillations**.

DHM[®] precisely samples up to 32'768 measurements along this waveform.




*cmUT membrane response to a **burst waveform composed of two periods of a sine wave at 5 MHz** (in blue). Vertical position versus time of the membrane at three monitoring area (in red, green and yellow) randomly selected over the measurement field of view (Courtesy to Philips).*

[Learn more about it](#)

Unique in the Frequency Domain: Mega-pixels Maps !


The Discrete Fourier Transform (DFT) of 3D time-sequences in the frequency domain **generalizes the conventional single point frequency analysis into mega-pixel vibration amplitude maps** calculated for each DFT bins/channels.



***cmUT membrane vibration amplitude map:** selection of three maps at 1.4, 3 and 4 MHz, among all those obtained by performing a DFT on each pixel of the 3D topography time-sequence response.*

[Learn more about it](#)

Characterize your MEMS in all environmental condition



In vacuum, at high or cryogenic temperatures, or in liquids, Lyncée proposes a full range of accessories and solutions ensuring:

- optimal optical imaging quality
- interferometric resolution
- fast and intuitive measurement
- friendly user software interface
- comprehensive data analysis

[Check out more](#)

Book a DHM[®] live-demo and help our forests!

Do you want to discover our product while avoiding unnecessary travels and interaction during the COVID-19 situation?

Until returning to normality, for any stand-alone DHM[®] purchased based on a [remote live-demo](#), Lyncée Tec commits to plant **200** trees through the non-profit organization [OneTreePlanted](#) and to provide you with a certificate.

[Contact us to book a live demo](#)



Lyncée Tec SA
Innovation Park, Bâtiment-A
CH-1015 Lausanne, Switzerland
info@lynceetec.com



This email was sent to
You received this email because you are registered with Lyncée Tec

[Unsubscribe here](#)

© 2020 Lyncée Tec