

DIGITAL HOLOGRAPHIC MICROSCOPY FOR COMPLEX PUPIL FUNCTION EVALUATION OF HIGH NA MICROSCOPE OBJECTIVES

Florian Charrière, Tristan Colomb and Christian Depeursinge
Ecole Polytechnique Fédérale de Lausanne (EPFL), Imaging and Applied Optics
Institute, CH-1015 Lausanne, Switzerland
E-mail: florian.charriere@a3.epfl.ch

KEY WORDS : Point spread function, pupil function, aberrations, digital holography




Figure 1: Axial (a) and radial (b) image comparisons in amplitude and phase: direct experimental APSF measurement (up); calculation based on pupil measurements (down).

proposed such an evaluation, based on pupil measurement with a Twyman-Green interferometer [3] but, for the first time to our knowledge, the accuracy of the pupil-based evaluation of the 3D APSF is directly compared (Fig. 1) with a quantitative measurement of the complex APSF function performed with a second original holographic system [4]. The validity of both techniques being definitely sealed through a comparison with a theoretical model [6]. Measurements on a 100x 1.3 NA MO are presented. The main advantage of the pupil evaluation regarding to the direct measurement techniques [4,6] resides in the single hologram required, having therefore no stability needs and a drastically reduced acquisition time.

Commonly, only the intensity point spread function is measured experimentally to describe an optical system, specifically a microscope objective (MO), neglecting the phase point spread function and losing therefore crucial information for aberrations quantification in the system. We present here a simple experimental method, which provides an accurate and reliable 3D evaluation of the complex amplitude point spread function (APSF) of a MO, based on a single holographic acquisition [1] of its pupil wavefront. The aberration function, extracted from the holographic pupil measurements, is developed in terms of Zernike polynomials and inserted in a scalar model of diffraction [2] to calculate the 3D distribution of the complex wavefront (amplitude and phase) around the focal point of the MO. Török already

[1] E. Cuche, P. Marquet and C. Depeursinge, "Simultaneous amplitude-contrast and quantitative phase-contrast microscopy by numerical reconstruction of Fresnel off-axis holograms", *Appl. Opt.* **38**, 6994-7001 (1999).

[2] M. Gu, *Advanced optical imaging theory*, Springer Verlag, 2000.

[3] P. Török and K. Fu-Jen, "Point-spread function reconstitution in high aperture lenses focusing ultra-short laser pulses", *Opt. Comm.* **213**, 97-102 (2002).

[4] A. Marian, F. Charrière, T. Colomb, F. Monfort, J. Kühn, E. Cuche, P. Marquet and C. Depeursinge, "On the complex 3D amplitude point spread function of lenses and microscope objectives. II Measurements by digital holography", submitted.

[5] S. F. Gibson and F. Lanni, *J. Opt. Soc. Am. A* **8**, 1601 (1991).

[6] R. Juskaitis and T. Wilson, "The measurement of the amplitude point spread function of microscope objective lenses", *J. Microsc.* **189**, 8-11 (1998).